



Pathology Center Memorandum: Change to Add-on Testing, 9/15/2020

Effective immediately, the 8 tests included in this document cannot be added on to tubes that have already had chemistry testing performed. This will ensure that we are complying with the suggestions of our reagent vendor, Roche Diagnostics, and will also ensure that we are not running a risk of 'carryover' in the chemistry analyzers, producing erroneous results.

We can still add-on/merge these tests to Extra Tubes and tubes that have not yet been analyzed.

BETA-HCG PREGNANCY

BETA-HCG TUMOR MARKER

CEA

HBSAG (DIALYSIS ONLY)

HIV - 1/2 Ag/Ab COMBO

LITHIUM

PCT

PSA

If you have questions, please contact our Client Support Services department at x44541. Thank you for your assistance in complying with this process.